



Patient Transfer Scale

for weighing time-critical or immobile patients safely, quickly and accurately.

Weighing critical patients is a crucial process and also a challenging one.

Weighing patients correctly is essential to ensuring accurate drug dosage and administering the correct treatment. For patients requiring urgent immediate medical attention, accurate weighing is not always easily achievable. It's for this reason, and through the determination of an acute nurse working for NHS Lanarkshire to find a reliable solution, that the Patient Transfer Scale was developed.

Key benefits include:

- ✓ One device for lateral transfer and weighing of patients
- ✓ Weigh immobile patients in seconds with less stress and more dignity for the patient.
- ▲ Accurate weight measurement when it's critical for safe medication management and patient monitoring



Instant Accurate Weighing when you need it most!

When knowing the weight of a patient is absolutely crucial to delivering the right care and treatment, the MDPTSM999 embeds confidence and accuracy in your processes.

Use the PTS throughout your hospital pathway to enhance outcomes particularly for; trauma patients in emergency departments where time is crucial, stroke patients to ensure a quick 'door to needle' time, plus Cardiac, ICU, HTU, and X-ray patients.

Time Saving Integrated Process

Our Patient Transfer Scale integrates a transfer board and patient weighing scale into one unit.

This means weighing patients becomes part of the already existing process of transferring a patient from a stretcher to bed, or bed to bed.

It's quick and easy, and saves you valuable time



Patient Transfer Scale

MDPTSM999



Easy to Read Display

The bright LCD is easy to read and can be programmed to automatically power-off after 30 or 60 minutes.

The backlight can easily be set to ON or OFF making the display even easier to read whatever lighting conditions you are in.

Lightweight and Portable

The MDPTSM999 has been made as lightweight as possible so it's easy to take it exactly where you need it and get an instant, accurate weight.

Grab handles are conveniently positioned down each side making it easier to move.

The scale is supplied with two hooks that conveniently fit into the grab handle space. This provides an easy method of wall mounting the scale when it's not in use.



Features Overview

- lightweight durable construction with low profile
 - easy to clean, wipe down surtace
 - 10 conveniently placed grab handles
 - weight hold function
 - bright easy to read LCD display
- display options include programmable auto-off & backlight
 - powered by rechargeable battery
- easily wall mount scale when not in use
- easy to use so requires minimal training
 - optional carry bag
 - WAND Registered


Freephone: 0800 800 379 enquiries@wedderburn.co.nz www.wedderburn.co.nz www.wedderburn.co.nz www.wedderburn.co.nz www.wedderburn.co.nz



Patient Transfer Scale

MDPTSM999

Features Overview



SPECIFICATIONS

| Legal for Trade Use | No |
|---|---|
| WAND Registered | Yes, Registration #190319-WAND-6S5XXQ |
| Australian Register of Therapeutic Goods (ARTG) | Yes, Registration # 136929 |
| Maximum Capacity | 250 kg |
| Graduation | 0.5 kg |
| Weighing Units | kg |
| Type of Display | 27.7 mm x 75 mm LCD |
| Environment Suitability | Dry |
| Function Keys | ON/OFF/ZERO, HOLD |
| Height Measurement | 100 cm to 175 cm (5 cm increments) |
| Operating Temperature | 5 to 35 °C |
| Power Source | Rechargeable battery pack (up to 50 hours continuous use) |
| Net Weight | 11.45 kg |
| Weighing Area (mm) | 1800 (W) x 510 (D) |
| Overall Dimensions (mm) | 1800 (W) x 700 (D) x 30 (H) |
| | |

OPTIONS

| Order Code | Description |
|-------------|-------------|
| MDPTSM999CB | Carry Bag |



Specifications may change without notice

Weighing Scales. Packaging Equipment. Labelling Equipment. Labels. Food Equipment. POS Solutions. Maintenance & Repairs.

Freephone: 0800 800 379